

A convenient carrying case, temperature cable, IBP extension cables, and your choice of three IBP adapters accompany the EasySim Remote Control and Educational Patient Simulator!

EasySim® Educational Patient Simulator Features:

- 12-Lead ECG
 - NSR
 - Tachycardia, Bradycardia, ST segment +/- 1 mm
 - ECG rates: +/- 10 BPM
 - PVC: Couplet, Run, Pause
 - VTach, VFib, Asystole
 - Bigeminy, Trigeminy, AFib
 - Atrial Flutter, 1st, 2nd, 3rd Degree Blocks
 - Pacer
 - Muscle and Baseline Noise
- Respiration
 - Adjustable Respiration Rates +/- 5 BPM
 - Apnea
 - Range: 0-120 BPM
- Invasive Blood Pressure: 4 Channels
 - Arterial, PA, CVP, LVP pressure Control
 - LV to AO Pullback Simulation (simulates transducer moving from LV to Ao and back)
 - CVP Swan Move Simulation (simulates transducer moving from CVP to RV to PA to Wedge and back)
 - Intra-aortic Balloon Pump (IABP)
 - User-controlled Variable Pressures
- Cardiac Output
 - 3, 5, and 7 liters/minute
 - Thermodilution CO Simulation
 - Injectate simulator not required
- Temperature YSI 700/400
 - 25, 37, 40 degrees C
- Unprecedented Battery Performance
 - No charging, no power cord, no changing batteries
 - Battery Runs Continuously for up to Five Years

Package Contents:

- (1) Multi-Parameter Arrhythmia Simulator
- (1) EasySim Remote
- (3) IBP Adapter Cables of your choice
- (3) IBP Extension Cables
- (1) Temperature Cable with stereo plug
- (1) Specially-modified nylon carrying case
- (1) 4-Year Warranty