

OX-2 OxSim Flex® Optical SpO2 Pulse Oximeter Simulator Features:

- Complete control of simulations!
- Saturations in 1% increments (10-100%)
- Pulse rates in 1 BPM increments (25-240)
- Create custom preset simulations for faster Alarm Testing
- Perfusion Index levels in 1% increments (10-100%)
- Factory Preset Simulations:
 - Saturation: 85%, Heart Rate: 80 bpm, Perfusion Index: approx. 2.0
 - Saturation: 95%, Heart Rate: 40 bpm, Perfusion Index: approx. 2.0
 - Saturation: 98%, Heart Rate: 80 bpm, Perfusion Index: approx. 2.0
 - Saturation: 98%, Heart Rate: 140 bpm, Perfusion Index: approx. 2.0
 - Saturation: 99%, Heart Rate: 80 bpm, Perfusion Index: approx. 2.0
- *Masimo Rainbow SET® Compatible*
- Automatic detection of SpO2 Manufacturer
- Color Display – OLED Technology
- Only 5 oz. (0.14 kg, excluding battery)
- Fast and Accurate
- Battery Operated – minimum 10 hours
- Power Save Mode
- Rugged: Polycarbonate housing. Drop tested 65 times from a height of 3 feet (1m) to a hard surface with no failures!
- Field-upgradable Software
- Four-year Limited Warranty

Package Contents:

- (1) OxSim® Flex® Optical SpO2 Pulse Oximeter Tester
- (1) Power Supply (6VDC)
- (1) AA battery
- (1) Certificate of Calibration
- (1) Carrying Case