

SC-3 SimCube® NIBP Simulator Features:

The SC-3 SimCube® NIBP Simulators are compact, affordable and easy to use. These hand-held simulators quickly and accurately determine whether or not a monitor's NIBP feature is working to spec. As cuff pressure changes, SimCube simulates real-life oscillometric pulses by varying both the size and shape of the wave. Get 150-250 NIBP simulations on compatible SimCube models using the optional Battery Boost with a set of 4 "AA" alkaline or NiMH batteries!

SC-3 SimCube® NIBP Simulator Features:

- Adult (120/80 mmHg), Neonatal (70/40 mmHg), and Hypertensive (190/120 mmHg) NIBP Simulations
- Digital Manometer
- Over Pressure (Peak Detect) Test mode measures the peak pressure produced by the patient monitor. The readout latches the highest pressure value sensed, in 0.5mmHg steps.
- Compatible with the Battery Boost Option

Package Contents:

- (1) SC-3 SimCube® NIBP Simulator
- (1) Power Supply (6VDC)
- (1) NIBP adapter
- (1) Certificate of Calibration

Optional Accessories: NIBP Adapters, Carrying Case, Cuff Jacket