

SC-4 SimCube® NIBP Simulator Features:

- Adult (120/80 mmHg), Neonatal (70/40 mmHg), and Hypertensive (190/120 mmHg) NIBP simulations
- Digital Manometer
- Over Pressure (Peak Detect) Test mode measures the peak pressure produced by the patient monitor. The readout latches the highest pressure value sensed, in 0.5mmHg steps.
- 12-lead ECG Simulation
- Respiration
- Pacer
- Heart Rate/Respiration sequence

Package Contents:

- (1) SC-4 SimCube® NIBP Simulator
- (1) Power Supply (6VDC)
- (1) Certificate of Calibration
- (1) NIBP adapter

Optional Accessories: NIBP Adapters, Carrying Case, Cuff Jacket, ECG Extender Module, Battery Boost