

SL-8 SimSlim® Multi-Parameter Patient Simulator Features:

- 12-lead ECG:
 - 8 Selectable Heart Rates: 30, 60, 90, 120, 140, 160, 200, 240 bpm
 - 10 Selectable Arrhythmias: NSR w/ PVCs, VTACH, VFIB, Asys, Bigeminy, Pacer, Trigeminy, ST+, 2nd Deg. Block, Square Wave
 - 8 Banana Plug Interface
- Respiration:
 - 8 Selectable Rates: 15, 30, 45, 60, 70, 80, 100, 120 bpm
- 4 IBP Channels:
 - 8 Dynamic: 120/80, 60/30, 30/10, 20/10, 125/15, 60/5, 8/2, 4/1 mmHg
 - 7 Static: 0, 25, 50, 100, 150, 200, 250 mmHg
- Cardiac Output Simulation:
 - 3, 5 and 7 l/m. No adapter cable required.
- YSI 700 and 400 Temp Simulations:
 - 25, 37, and 40 degrees Celsius
- Battery Life:
 - Total life depends on usage: 5-10 years, all day, every day. No charging or changing required.

Package Contents:

- (1) SL-8 SimSlim Simulator
- (1) IBP adapter of customer choice
- (1) IBP extension cable, 6ft.
- (1) Temperature adapter cable
- (1) Certificate of Calibration

Optional Accessories: Additional IBP Adapters