



Safe-T Sim® Mobilize
Safety Analyzer



FlowTrax® Mobilize
IV Pump Analyzer/Press. Meter



SimCube® Mobilize
ECG, Resp, IBP,
NIBP Simulator



OxSim Flex® Mobilize
SpO₂ Simulator



TOOLS OF THE TRADE

Pronk Technologies

Mobilize — Wireless, Automated Test Equipment

Mobilize is a complete wireless solution that connects Safe-T Sim® Safety Analyzer, SimCube® and OxSim Flex® Vital Signs Simulators and FlowTrax® IV Pump Analyzer to your smart device. This provides complete flexibility to control all these products from one user interface, set up your own test protocols or run your CMMS procedures right from your Apple iOS device.

Key features include: instantaneous wireless synchronization to Mobilize products, including all simulations, as well as full control of all IV pump and electrical safety analyzer testing functions via the Mobilize Apps. Capture all tests performed; measurements are recorded automatically while you work; or simply tap Pass/Fail. Add images of device under test into a complete detailed electronic report. App integrated Barcode Scanner and user-configurable wireless PM Sticker printing. Use app autonomously or integrate to a CMMS. Stores up to 1000 detailed test records with *no* monthly/annual license fee required. ⚙️

For more information, visit www.pronktech.com/mobilize.